

Anytime Align






TABLE OF CONTENTS

Our Mission	01
Три	03
Why Anytime Align	04
Product Specifications	08
Quality Assurance	09
Case Types	10
Aligner Solutions	11
Pricing	13
Packaging Designs	14
Ordering Process	16
Treatment	17
Planning Process	10
Support And Training	18
Our Products	19
Case Studies	21
Our Promise	25







TPU is renowned for its unique combination of flexibility, durability, and biocompatibility, making it the ideal choice for orthodontic treatment."



ANYTIME BUTTONS TECHNOLOGY

PLANNING SOFTWARE

Access to millions of case data points.

Al-driven treatment simulations.

Proprietary anchor points for enhanced treatment control. Seamless integration with elastics.



EXCLUSIVE TECHNICIAN ASSIGNMENT

Personalized service for your practice. Consistent point of contact for all cases.



TPU Material:

Precision, Comfort, and Safety in Every Aligner

At the heart of our aligners lies Thermoplastic Polyurethane (TPU), a cutting-edge material engineered to deliver exceptional performance, comfort, and reliability.

TPU is renowned for its unique combination of flexibility, durability, and biocompatibility, making it the ideal choice for orthodontic treatment.

Why TPU?

Superior Comfort:

TPU is soft yet resilient, ensuring a snug fit that minimizes discomfort while providing consistent pressure to move teeth effectively.

Durability:


Its high resistance to wear and tear ensures that each aligner maintains its shape and effectiveness throughout the treatment period.

Precision:

TPU's flexibility allows for precise tooth movement, enabling us to achieve optimal treatment outcomes tailored to each patient's unique needs.

Transparency:

The material is crystal clear, offering a discreet solution for those seeking an invisible orthodontic treatment.



WHY ANYTIME ALIGN

STAND ABOVE THE REST





More Transparent

Designed for near-invisibility, Anytime Align offers a discreet and natural look. Plus, they're more comfortable and resist stains better than other leading aligner brands.






Enhanced for a Better Fit

Crafted with precision-polished and contoured edges, Anytime Align ensure a smoother and more comfortable fit for an improved treatment experience.





Engineered for Stain Resistance.

Our advanced material is designed to resist discoloration, keeping your aligners looking clearer and fresher throughout your treatment.



BEFORE AFTER PATIENT ALIGNER











AFTER



BEFORE



Safety You Can Trust

Our TPU material is FDA and CE certified, meeting the highest standards of safety and quality. These certifications guarantee that our aligners are biocompatible, free from harmful substances, and suitable for prolonged use in the mouth.

Commitment to Excellence

By choosing TPU, we ensure that every aligner is not only effective but also safe and comfortable for our patients. This commitment to quality is why our aligners stand out in the market, providing a reliable and innovative solution for achieving the perfect smile.

With TPU, you're not just getting an aligner—you're getting a product designed with precision, care, and your well-being in mind

Modular Treatment Approach

- 10-Aligner Shipments: Aligners are sent in batches of 10 at a time, allowing real-time adjustments.
- Dynamic Material & Treatment Adaptation: Unlike traditional providers, we optimize both material and treatment as the case progresses.
- Before the 10th aligner is placed, a new intraoral scan is conducted, allowing the next set to be prepared and delivered 5 days in advance, ensuring seamless transitions, more accurate tooth movements, fewer delays, and higher treatment success rates by continuously tailoring the plan to the patient's real-time progress.
- Ethnicity-Based Treatment Planning
- Beyond Big Databases: More data doesn't always mean better outcomes due to the law of diminishing returns in Al. Instead, ethnicity-specific anatomical variations play a key role in achieving optimal treatment results.



Our Guarantee

Fast Turnaround: Treatment plans delivered within 2 weeks, or receive a discount on the case




Why It Matters:

- Bone Structure Variations Jaw shape, arch width, and bone density differ across ethnic groups, requiring customized treatment strategies.
- Soft Tissue & Lip Support Ethnicity influences lip thickness and soft tissue response, affecting esthetic outcomes.
- Tooth Size & Crowding Asian populations are more prone to dental crowding, necessitating precise space management and movement strategies.
- Cultural Aesthetic Preferences Aligners are designed to match unique esthetic expectations for a natural, personalized smile.

How We Optimize:

- Al-Trained on Ethnic Data for improved movement predictions.
- Adaptive Treatment Adjustments to enhance efficiency.
- Specialized Crowding Solutions for better space management and alignment.
- By prioritizing quality over quantity in data, we ensure more precise, efficient, and esthetically superior treatment—setting us apart from traditional aligner providers.





CASE TYPES

SMILES FOR EVERYONE!

Clear aligners are a simple and effective way to correct a variety of dental issues. Their nearly invisible design makes them a popular choice, as they blend in naturally and encourage consistent wear for better results.



Satisfaction guarantee

Comprehensive warranty policy






Our subscription-based model makes clear aligner treatment more accessible and affordable. By streamlining the process, we help dental professionals deliver high-quality care with ease, ensuring patients achieve straighter smiles effortlessly.

Industry-Leading Expertise:

Our team consists of industry veterans with proven experience from top aligner companies like Byte, Invisalign, and AngelAlign, backed by prestigious certifications such as GDC and global registrations.



Revolutionary **End-to-End Solutions:**

We offer a complete suite of services, from advanced treatment planning and production overseen by a global network of orthodontists to private-label aligners with custom branding, packaging, and clinic management software.

Beyond Aligners:

Our services extend to free design of veneers, crowns, and dentures, website development, social media management, and catalog design—all included to help clinics reduce chair time, minimize refinements, and cut costs.

Proven Results:

Multiple clients have experienced a 25% decrease in costs with our 1-year package and a 15% reduction with our 6-month package, thanks to our full range of services and cutting-edge technology.

Empowering Clinics:

We help practices streamline operations, enhance patient satisfaction, and significantly boost profitability, enabling them to thrive in a competitive landscape.

Anytime Align Advantage:

Revolutionize your practice with our subscription-based, white-label clear aligner service. We provide cutting-edge technology, 24/7 support, and dedicated technicians, empowering you to deliver top-tier orthodontic care under your own brand while growing your practice.

Orthodontics, Reimagined: With our world-class expertise, innovative tools, and holistic approach, we're not just delivering aligners—we're reshaping the future of orthodontic care.

> Contact us today to learn how we can transform your practice and help you achieve similar success!

PRICING

Our Innovative Based Pricing

For a 6-month aligner treatment, each set costs MYR 135, totaling MYR 2,700 for 20 sets.

With our subscription plan, you get free cases and complimentary lab designing for 3 months. The subscription rate is MYR 145 per set, bringing the total to MYR 2,900.

For a 12-month aligner treatment, each set costs MYR 125, totaling MYR 2,500 for 20 sets.

A single case without a subscription costs MYR 165 per set, making it MYR 3,300 for 20 sets.

This plan offers added value, making it a cost-effective choice for practitioners.

165 MYR Single case **125 MYR** 12 months

135 **MYR** 6 months

145 MYR 3 months

BRONZE RM 1,999 / month 3 Months Lab services 4per month Free Delivery Premium Packaging (custom logo only) Free Refinement **Business Hour Customer Services** Free Retainers

RM 5,999

SILVER RM 1,499 / month 6 Months Lab services oper month Free Delivery Premium Packaging (custom logo only) Free Refinement/Retainer 24/7 Customer Service 1 free case Worth RM 5,000 1 case 40% Off RM 8,999

GOLD **RM 1,333** / month 12 Months Lab services 8per month Free Delivery Premium Packaging Custom logo with options for size, color scheme, additional personalization (pattern, interior lining color, packaging box quality), engraving, and embossing. Free Refinement 24/7 Customer Service Two Complimentary Cases First Single/ Double Case Purchase (40% Off) RM 15,999

Free clinic management softwaree basic for bronze medium for sliver and advanced for gold

ELEGANT PACKAGING DESIGNS



Soft and Warm Vanilla Dream



Elegant and Classic **Midnight Velvet**



Modern and Sleek **Tuxedo Chic**



Ocean Depths **Calm and Serene**



Sweet and Cheerful **Blossom Pink**



Fresh and Clean **Snowy Whisper**





All words on the pouch is customisable the pouch is made out of aluminium and is vacum sealed



Color Option




Picture of inside the box







OUR EFFORTLESS ORDERING PROCESS



TREATMENT PLANNING PROCESS



SUPPORT AND TRAINING Initial onboarding and software training · Ongoing education and webinars · 24/7 customer support · Dedicated account manager Private Label Clear Aligner Solutions

Our products

Precision, Comfort, and Safety in Every Aligner



ALIGNERS WITH ANYTIME ANCHORS



KIDS ALIGNERS



SNORE GUARDS



ANYTIME MYO



BEFORE AFTER PATIENT ALIGNER











AFTER



BEFORE



CASE STUDIES


AFTER WEARING 35 PAIRS OF ALIGNERS


BEFORE AFTER



AFTER WEARING 45 PAIRS OF ALIGNERS

BEFORE





ALIGNING ARCHES &

CORRECTING DEEP BITE IN 40 STEPS

This case presents malaligned upper and lower arches, bimaxillary proclination, and a deep bite. Correction is being achieved through distalization of the maxillary posterior segment, followed by retraction of the maxillary anteriors. Deep bite is addressed with bite turbos placed on the lingual surface of the maxillary anteriors and bite ramps added to the occlusal surface of the mandibular molars.

Additionally, EBTs have been designed to align the mandibular left second premolar. The treatment plan consists of 40 steps, focusing on improving overjet, overbite, and overall arch alignment.


BEFORE















OUR CASE DESIGN

SOFTWARE AND MANAGEMENT SOFTWARE



CASE STUDIES

AFTER WEARING 43 PAIRS OF ALIGNERS

BEFORE









AFTER







GETTING STARTED



OUR PROMISE

- Transparent, subscription based pricing
- Volume discounts available

No hidden fees

Flexible payment options



CONTACT US

Phone Number: +601123453772 Email: info@anytimealign.com Website: www.anytimealign.com



